



US006752810B1

(12) **United States Patent**
Gao et al.

(10) **Patent No.:** **US 6,752,810 B1**
(45) **Date of Patent:** **Jun. 22, 2004**

(54) **INSTRUMENT AND METHOD FOR PULLING AND TWISTING A TIE ONTO TWO SEPARATED ITEMS**

(75) Inventors: **Hua Gao**, Fox Point, WI (US); **James A. Rinner**, Racine, WI (US)

(73) Assignee: **Pilling Weck Incorporated**, Horsham, PA (US)

(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 106 days.

(21) Appl. No.: **10/158,470**

(22) Filed: **May 31, 2002**

(51) **Int. Cl.**⁷ **A61B 17/56**

(52) **U.S. Cl.** **606/103; 606/74**

(58) **Field of Search** 606/103, 74, 148; 140/93 R, 93.2, 117, 118, 119

(56) **References Cited**

U.S. PATENT DOCUMENTS

1,120,575 A	12/1914	Wertz	
1,209,434 A	12/1916	Hayden	
1,304,620 A	5/1919	Steinkoenig	
1,365,649 A	1/1921	Bates	
1,463,869 A	8/1923	Campbell	
2,049,361 A	7/1936	Ericsson	
2,279,068 A	4/1942	Siebrandt	
2,291,413 A	7/1942	Siebrandt	
2,455,609 A	12/1948	Schieb	
2,657,718 A	11/1953	Greathouse	
3,273,605 A	9/1966	Ferrara	
3,507,270 A *	4/1970	Ferrier	600/481
3,759,302 A	9/1973	Attenborough	
3,865,155 A	2/1975	Spath	
4,512,346 A	4/1985	Lemole	
4,527,554 A	7/1985	Klein	
4,587,963 A	5/1986	Leibinger	
4,606,335 A	8/1986	Wedeen	
4,813,416 A	3/1989	Pollak	
4,880,038 A	11/1989	Meinershagen	
4,935,027 A *	6/1990	Yoon	606/146
4,966,600 A	10/1990	Songer et al.	
5,004,020 A	4/1991	Meinershagen	

5,116,340 A	5/1992	Songer	
5,127,448 A	7/1992	Hillestad	
5,275,383 A	1/1994	Wick	
5,312,410 A	5/1994	Miller	
5,417,698 A	5/1995	Green	
5,447,512 A	9/1995	Wilson	
5,449,361 A	9/1995	Preissman	
5,501,688 A	3/1996	Whiteside	
5,540,698 A	7/1996	Preissman	
5,601,572 A	2/1997	Middleman	
5,741,279 A *	4/1998	Gordon et al.	606/144
5,752,551 A	5/1998	Trueblood	
5,772,663 A	6/1998	Whiteside	
5,810,832 A	9/1998	Blasingame	
5,830,234 A	11/1998	Wojciechowicz	
5,849,012 A	12/1998	Abboudi	
5,851,209 A	12/1998	Kummer	
5,935,133 A *	8/1999	Wagner et al.	606/103
5,968,077 A	10/1999	Wojciechowicz	
6,045,551 A	4/2000	Bonutti	
6,045,561 A	4/2000	Marshall	
6,185,858 B1	2/2001	Choron	
6,217,580 B1	4/2001	Levin	
6,287,307 B1	9/2001	Abboudi	
6,332,889 B1	12/2001	Sancoff	
D476,084 S *	6/2003	Gao et al.	D24/133
2003/0028202 A1 *	2/2003	Sancoff et al.	606/144

* cited by examiner

Primary Examiner—Michael J. Milano
Assistant Examiner—Bradford C. Pantuck
(74) *Attorney, Agent, or Firm*—Arthur J. Hansmann

(57) **ABSTRACT**

An instrument and method for applying a suture of stiff but bendable material, such as wire, to broken bone for closing and holding the break. Also, the instrument and method are for pulling and holding two items together with a tie. The instrument has two parts movable relative to each other for pulling the items or bone together and for twisting the suture onto the bone to hold it together. The suture is strung onto the instrument on the exterior thereof without having the suture extend through the instrument interior, so the surgeon can guide the suture and accurately string it. A crank is applied to twist the suture after the suture has been initially secured to the bone to hold the bone together.

22 Claims, 5 Drawing Sheets

